



WOUND CARE
BECAUSE YOU CARE



FOAM NON-ADHESIVE

WOUND CARE
BECAUSE YOU CARE


FEATURES

Excellent absorption of exudate // Promotes healing through a moist wound environment // Reduces the risk of maceration // Low friction backing // Low adherent wound contact layer // Soft and conformable // Ideal under compression //

ACTIVHEAL® FOAM NON-ADHESIVE CREATES A MOIST WOUND HEALING ENVIRONMENT


ActivHeal® Foam Non-Adhesive is a versatile dressing developed to offer protection and absorbency on moderate to heavily exuding wounds.

Its three-layer construction offers a low friction backing, a soft highly absorbent foam and a perforated wound contact layer. This non-adhesive dressing is ideal when reducing trauma is a key part of the wound dressing regime. It has soft and flexible properties mean patients can wear it underneath clothing with comfort and confidence. The ActivHeal® Foam Non-Adhesive range offers a variety of shapes and sizes including; the ActivHeal® Foam Heel dressing, a two layer, soft pre-moulded dressing, specifically designed to apply easily to the heel and other bony anatomical areas. The ActivHeal® Foam Tracheostomy two layer dressing features a versatile fenestration, allowing it to fit neatly around a tracheostomy tube, intubation tube, cannula insertion, external bone fixators and appropriately sized wound drainage tubes.



PERFORMANCE

Total Fluid Handling (g/cm²)¹



SIZES AND CODES

SIZE (CM)	DRESSINGS PER CARTON	PRODUCT CODE	NHS SUPPLY CHAIN CODE	DT PIP CODE
5x5	10	10009113	ELA214	314-6073
7.5x7.5	10	10009114	ELA245	NA
10x10	10	10009115	ELA216	304-3791
10x20	10	10011265	ELA246	317-0503
20x20	10	10009117	ELA218	314-6081
FOAM HEEL				
18x12	5	10010227	ELM160	356-1131
FOAM TRACHEOSTOMY				
10x10	10	10009118	ELA335	382-2483

CASE STUDY

Patient A was admitted to the unit with a fresh, sizeable skin tear to her left shin. The skin on her right leg was in danger of further breakdown. Her past medical history included left hemiparesis due to cerebral infarct, osteoarthritis, hypertension, shingles and low mood.



- Day 1 - The wound was assessed and the skin tear measured 7cm at the widest point. The dermis and the epidermis had been pulled apart. Distally the epidermal flap was lost, leaving a small area of dermal tissue exposed. The wound had no areas of necrosis or slough, moderate levels of exudate, and no sign of infection. Additionally the flap was also viable.
- Day 11 - The wound continued to show improvement following the use of ActivHeal® Non-Adhesive Foam. Proximal wound margins were now apposed. No signs of maceration or breakdown of surrounding skin. The wound area had reduced in size and contained both granulation tissue and epithelial tissue.
- Day 31 - The wound showed a good level of healing and, as a result, the patient was discharged from the unit back into the community. The District Nurses were requested to moisturise both legs, then continue with ActivHeal® Non-Adhesive Foam, held in place with a tubular bandage to protect the delicate, newly healed tissue and prevent tissue breakdown.

INDICATIONS

Activ-Heal® Foam Non-Adhesive is indicated for moderate to heavily exuding wounds:

- Pressure ulcers
- Venous Leg ulcers
- Diabetic ulcers
- Post operative surgical wounds
- Superficial and partial thickness burns
- Lacerations
- Abrasions
- Graft wounds and Donor sites
- Cavity wounds as a secondary dressing

Moderate > Heavy Exudate



Wound Types



Advanced Medical Solutions Ltd

Premier Park, 33 Road One, Winsford Industrial Estate, Winsford, Cheshire CW7 3RT
T 08444 125 755 **E** activheal@admedsol.com **W** www.activheal.com



References 1. P3422R

ActivHeal, its logo and the Advanced Medical Solutions logos are registered trademarks of Advanced Medical Solutions Ltd
 © Advanced Medical Solutions 2020 // Review date May 2020 // Document code AMS/AH-004-02.